


ICT+ 30ML STERILE 70% IPA POUCH

SATURIX

WFI Quality Wate

Precision Dosing for Isolator Cleaning

Consistent and precise saturation of every ICT+ pad to improve compliance and minimize chemical waste.



- 70% USP Grade Isopropyl Alcohol
- 30% Water for Injection (WFI) quality
- Alcohol blend filtered to .1 microns
- Pyrogen free
- Gamma irradiated and validated to a sterility level of 10⁻⁶

 Cleanroom Level Recommendation*

 ISO
 3
 4
 5
 6
 7
 8

Item	Size	Quantity	
IPA30S	30ml	4 pouches / bag, 15 bags / case	

* Cleanroom class recommendations are based on typical industry-wide applications and general cleanliness classifications. Please consult our experienced FG Clean Wipes sales epresentatives for specific cleanroom selections for critical and/or specialized applications based on a wider selection of criteria. To learn more about our products, visit us online at Saturix.com or contact us at +1 (413) 598-8344 or info@saturix.com

